



7018 Owensmouth Ave. Suite 103
 Canoga Park, CA, 91303
 Phone: 818-710-1281
 Fax: 818-936-0121


Email: Info@immunospec.com
www.immunospec.com



Trichinella IgG



Catalog No: E21-284



FOR INVITRO DIAGNOSTIC USE ONLY

Intended Use

For the qualitative screening of serum IgG antibodies to *Trichinella* using an Enzyme Linked Immunoabsorbant Assay (ELISA) technique.

Summary

Immunospec Trichinosis, the infection caused by the nematode *Trichinella spiralis*, is acquired by ingestion of raw or undercooked meats (primarily pork). Although the nematode may be found in a wide variety of animals worldwide, the domestic pig is the primary source of infection in developed nations.

Serology has also been an important tool in the diagnosis of trichinosis for several decades. Various methodologies; such as ELISA, latex agglutination (LA), indirect hemagglutination (IHA) and bentonite flocculation (BFT) have been used. Although various classes of antibodies have been detected, no single class has shown superior diagnostic ability over the others.

BFT has been the method of choice for serology but suffers from nonspecific reactions, some lack of sensitivity (measurable antibodies often do not appear until 3 to 4 weeks after infection) and difficulty in performing the test. Recently, an excretory-secretory (ES) antigen has been purified from the larvae of infected pigs. This antigen has a high degree of specificity for *T. spiralis* and has been used in several large scale studies.

Principle of Procedure

The micro test wells are coated with *Trichinella* antigen. During the first incubation with the diluted patients' sera, any antibodies which are reactive with the antigen will bind to the coated wells. After washing to remove the rest of the sample, the Enzyme Conjugate is added. If antibodies have been bound to the wells, the Enzyme Conjugate will then bind to these antibodies. After another series of washes, a chromogen (tetramethylbenzidine or TMB) is added. If the Enzyme Conjugate is present, the peroxidase will catalyze a reaction that consumes the peroxide and turns the chromogen from clear to blue. Addition of the Stop Solution ends the reaction

and turns the blue color to a bright yellow color. The reaction may then be read visually or with an ELISA reader.

Reagents

Item	Description
Test Strips	Microwells containing <i>Trichinella</i> antigens – 96 test wells in a test strip holder.
Enzyme Conjugate	One (1) bottle containing 11 ml of Protein A conjugated to peroxidase.
Positive Control	One (1) vial containing 2 ml of diluted positive rabbit serum.
Negative Control	One (1) vial containing 2 ml of diluted negative human serum.
Chromogen	One (1) bottle containing 11 ml of the chromogen tetramethyl benzidine (TMB).
Wash Concentrate (20X)	One (1) bottle containing 25 ml of concentrated buffer and surfactant.
Dilution Buffer	Two (2) bottles containing 30 ml of buffered protein solution.
Stop Solution	One (1) bottle containing 11 ml of 1 M phosphoric acid.

Precautions

Do not use solutions if they precipitate or become cloudy. Wash concentrate may show crystallization upon storage at 2 – 8 °C. Crystallization will disappear after dilution to working strength.

Do not use serum that may have supported microbial growth, or is cloudy due to high lipid content. Samples high in lipids should be clarified before use.

Treat all sera as if capable of being infectious. Negative control has been tested and found negative for Hepatitis B surface antigen and for the antibody to HIV by required test methods. This product should be used under appropriate safety conditions that would be used for any potentially infectious agent.

Do not add azides to the samples or any of the reagents.

Storage Conditions

Reagents, strips and bottled components:

Store between 2 – 8 °C.

Squeeze bottle containing diluted wash buffer may be stored at room temperature.

Preparation

Wash Buffer - Remove cap and add contents of bottle to 475 ml of reagent grade water. Place diluted wash buffer into a squeeze bottle with a narrow tip opening.

Note: Washings consist of filling to the top of each well, shaking out the contents and refilling.

Avoid generating bubbles in the wells during the washing steps.

Collection and Preparation Of Serum

Serological specimen should be collected under aseptic conditions. Coagulate blood and remove serum. Hemolysis is avoided through prompt separation of the serum from the clot.

Serum should be stored at 2-8 ° if it is to be analyzed within a few days. Serum must be held for 3 to 6 months by storage at -20 °C. Lipemic and strongly hemolytic serum should be avoided. Do not heat inactivate serum and avoid repeated freezing and thawing of samples.
 Test samples: Make a 1:64 dilution of patient's sera using the dilution buffer (e.g. 5 µl sera and 315 µl dilution buffer).

Procedure

Materials Provided

Trichinella Serology Microwell ELISA Kit

Materials Required But Not Provided

Pipettes
 Squeeze bottle for washing strips (narrow tip is recommended)
 Reagent grade water and graduated cylinder
 Tubes for sample dilution
 Absorbent paper

Suggested Materials

ELISA plate reader with a 450 nm and a 650 to 620 nm filter (optional if results are read visually)

Performance of Test

- Break off number of wells needed (two for controls plus number of samples) and place in strip holder.
- Add 100 µl (or two drops) of the negative control to well #1, 100 µl of the positive control to well #2 and 100 µl of the diluted (1:64) test samples to the remaining wells.
 Note: Negative and positive controls are supplied prediluted. Do not dilute further.
- Incubate at room temperature (15 to 25 °C) for 10 minutes.
- Shake out contents and wash 3 times with the diluted wash buffer.
- Add 2 drops of Enzyme Conjugate to each well.
- Incubate at room temperature for 5 minutes.
- Shake out contents and wash 3 times with wash buffer. Slap wells against paper towels to remove all of the wash buffer.
- Add 2 drops of the Chromogen to every well.
- Incubate at room temperature for 5 minutes.
- Add 2 drops of the Stop Solution and mix by tapping strip holder.

Reading of Results

Visually: Look at each well against a white background (e.g. paper towel) and record as clear or +, ++ or +++ reaction.
 ELISA Reader: Zero reader on air. Set for bichromatic readings at 450/650-620 nm.

Test Limitations

Serologic results are an aid in diagnosis but cannot be used as the sole method of diagnosis.

Quality Control

The use of controls allows validation of kit stability. The kit should not be used if any of the controls are out of range. Expected values for the controls are:

Negative - 0.0 to 0.3 OD units
Positive - 0.5 OD units and above

Troubleshooting

Negative control has excessive color after development.

Reason: inadequate washings.

Correction: wash more vigorously. Remove excessive liquid from the wells by tapping against an absorbent towel. Do not allow test wells to dry out.

Interpretation of Results - ELISA Reader

Zero ELISA reader on air. Read all wells at 450/650 to 620 nm.

Positive - Absorbance reading equal to or greater than 0.3 OD units.

Negative - Absorbance reading less than 0.3 OD units.

A negative OD reading indicates that the patient has no detectable level of antibodies. This may be due to lack of infection or poor immune response by the patient.

Interpretation of Results -Visual

Compare results to the controls. A sample should be interpreted as positive if the degree of color development is obvious and significant.

Expected Results

The number of individuals showing positive results can vary significantly between populations and geographic regions. If possible, each laboratory should establish an expected range for its patient population.

Performance Data

Study #1- Canadian Reference Center

Compared IMMUNOSPEC ELISA to another commercial ELISA. Found concordance of 85.4% (n=82).

Study #2 - CDC&P

		CDC&P	
		+	-
IMMUNOSPEC	+	43	1
	-	2	15

Sensitivity (biopsy positive) of 94.4% (17/18)

Sensitivity (outbreak symptomatic) of 96.3% (26/27)

Specificity - Normals 93.8% (15/16)

Ascaris 100% (6/6)
 Hookworm 83.3% (5/6)
 Strongyloides 83.3% (5/6)
 Toxocara 66.6% (4/6)
 Trichuris 83.3% (5/6)

References

- Despommier, D.: Trichinellosis. Immunodiagnosis of Parasitic Diseases, Vol. 1, *Helminthic Diseases*. Ed. Walls and Schantz. Academic Press, 1986. pp. 163-181.
- Krogstad D., Visvesvara G., Walls R., Smith J.: Tissue Helminths. *Manual of Clinical Microbiology*, 4th Ed.. American Society for Microbiology, Washington, DC. 1985, pp. 654.
- Kagan, I.: Serodiagnosis of Parasitic Diseases. *Manual of Clinical Laboratory Immunology*. 3rd Ed.. American Society for Microbiology, Washington, DC. 1986, p. 474-477.
- Petri, W., et.al.: Common-Source Outbreak of Trichinosis Associated With Eating Raw Home-Butchered Pork. *So Medical J*, August 1988, pp. 1056-1058.

5. Oliver, D., et.al.: Enzyme Linked Immunoassay For Detection of Trichinosis In Humans. *Proc 7th Int Conf on Trich*, Oct. 1988.
6. Ivanoska, D., et.al.: Comparative Efficacy of Antigen and Antibody Detection Tests for Human Trichinellosis. *J Parasit*, Vol. 75 #1, 1989, pp. 38-41.



European Authorized Representative:

CEpartner4U , Esdoornlaan 13, 3951DB Maarn
. The Netherlands. Tel.: +31 (0)6.516.536.26



Manufacturer:

IMMUNOSPEC CORPORATION
7018 Owensmouth Ave. Suite 103
Canoga Park, CA, 91303
Phone: 818-710-1281
Fax: 818-936-0121

PI E21-284
Revision B
Revision Date: 12/29/2010