



7018 Owensmouth Ave. Suite 103
Canoga Park, CA 91303
Tell : (818)710-1281
Fax : (818)936-0121
Email: Info@immunospec.com
www.immunospec.com



Hantavirus IgG/IgM Combo Test

Cat.# R13-524

A Rapid qualitative Immunochromatographic test for the simultaneous detection of IgG and IgM Antibodies to Hantavirus in human Whole Blood, Serum or plasma

INTRODUCTION

Immunospec Hantavirus Test is a rapid immunochromatographic assay for the simultaneous detection of IgG and IgM antibodies to Hantavirus in human whole blood, serum or plasma. The assay is used as a screening test for Hantavirus infection and as an aid for differential diagnosis of primary and secondary infections in conjunction with other criteria.

SUMMARY

Hantaviruses are rodent-borne pathogens and are normally transmitted to humans via aerosols generated from feces, urine, and saliva of infected rodents. Hantaviruses are the major cause of human Hemorrhagic Fever with Renal Syndrome (HFRS). HFRS is a group of clinically similar illnesses, including diseases such as Korean hemorrhagic fever, epidemic hemorrhagic fever and nephropathic epidemics. Hantaviruses belong to family *Bunyaviridae*. The subtypes that cause HFRS include Hantaan, Dobrava-Belgrade, Seoul and Puumala. HFRS is found throughout the world. Symptoms of HFRS usually develop within 1 to 2 weeks after exposure to infectious material, but in rare cases, they may take up to 8 weeks to develop. Initial symptoms begin suddenly and include intense headaches, back and abdominal pain, fever, chills, nausea, and blurred vision. Individuals may have flushing of the face, inflammation or redness of the eyes, or a rash. Later symptoms can include low blood pressure, acute shock, vascular leakage, and acute kidney failure which can cause severe fluid overload. The severity of the disease varies depending upon the virus causing the infection. Hantaan and Dobrava virus usually cause severe symptoms, while Seoul and Puumala virus infections are usually more moderate. Complete recovery can take weeks or months. Depending upon which virus is causing the HFRS, death occurs in less than 1% to as many as 15% of patients. The early detection of hantavirus infection can help patient to be treated promptly to avoid the severe result. Specific antibody responses to Hantavirus virus enable serodiagnosis and differentiation between primary and secondary dengue infections.

Immunospec Hantavirus Test is a new generation rapid Immuno-chromatographic test using recombinant viral antigens to detect specific antibody response.

TEST PRINCIPLE

Immunospec Hantavirus Test utilizes the principle of Immunochromatography. Anti-human IgM and human IgG antibodies are immobilized on the nitrocellulose membrane respectively, as two individual test lines (IgM line and IgG line) in the test window of the test device. The IgG line in the test window is closer to the sample well and followed by IgM line. As the test sample flows through the membrane assembly within the test device, the colored-Hantavirus specific recombinant antigen-colloidal gold conjugate complexes with specific antibodies (IgM or IgG) of Hantavirus virus, if present in the sample. This complex moves further on the membrane to the test region where it is immobilized by the anti-human IgM and/or human IgG binding proteins coated on the membrane leading to formation of a colored band, which confirms a positive test results. Absence of this colored band in the test window indicates a negative test result. A built-in control line will always appear in the test window when the test has performed properly, regardless of the presence or absence of anti-Hantavirus virus antibodies in the specimen.

REAGENTS AND MATERIALS SUPPLIED

Each kit contains:

1. Immunospec Hantavirus Test card in foil pouch
2. Sample buffer
3. Product insert

MATERIALS NOT PROVIDED

1. Specimen collection container
2. Timer

STORAGE AND STABILITY

The sealed pouches in the test kit may be stored between 2-30°C for the duration of the shelf life as indicated on the pouch.

PRECAUTIONS

1. This kit is for *IN VITRO* diagnostic use only.
2. This kit is for **PROFESSIONAL** use only.
3. Read the instructions carefully before performing the test.
4. This product does not contain any human source materials.
5. Do not use kit contents after the expiration date.
6. Handle all specimens as potentially infectious.
7. Follow standard Lab procedure and biosafety guidelines for handling and disposal of potentially infective material. When the assay procedure is complete, dispose specimens after autoclaving them at 121°C for at least 20 min. Alternatively, they can be treated with 0.5% Sodium Hypochlorite for 1-2 hours before disposal.
8. Do not pipette reagent by mouth and no smoking or eating while performing assays.
9. Wear gloves during the whole procedure.

SPECIMEN COLLECTION AND PREPARATION

1. No prior special preparation of the patient is required before sample collection by approved techniques.
2. Fresh serum / plasma is preferable. Serum / plasma may be stored at 2-8°C up to 3 days in case of delay in testing. For long-term storage, freeze the specimen at -20°C for 3 months or -70°C for longer periods.
3. The test works best on fresh whole blood samples. If testing cannot be done immediately, Blood samples collected with a suitable anticoagulant such as EDTA or Heparin or Oxalate should be stored at 2-8°C up to 3 days. Blood samples should not be frozen.
4. Repeated freezing and thawing of the specimen should be avoided.
5. Do not use haemolysed, clotted, contaminated, lipamic and viscous/turbid specimen.
6. Specimen containing precipitates or particulate matter must be centrifuged and the clear supernatant only used for testing.
7. Do not heat inactivate the sample.
8. Shipment of samples should comply with local regulations for transport of etiologic agents.

PROCEDURE

1. Bring the kit components to room temperature before testing.
2. Open the pouch and remove the Card. Once opened, the test card must be used immediately.
3. Label the test card with patients identity.
4. Apply 5 µL of serum, plasma or whole blood to the "S1" area indicated by the **arrow mark**.
5. Add 2 to 3 drops of sample buffer to sample well marked as "S".
6. At the end of 10 minutes read the results. A strong positive sample may show result earlier.

INTERPRETATION OF RESULTS

1. Negative
Only control line appears.
2. IgM Positive
Both control line and the second test line (the higher test line) appear. It indicates the possibility of primary infection.
3. IgM and IgG Positive
Control line and both test lines appear. It indicates the possibility of acute secondary infection.
4. IgG Positive
Both control line and the second test line (the lower test line which is closer to the sample well) appear. It indicates the possibility of the secondary infection or past infection.

5. Invalid Result

If after 10 minutes no line is visible within the test or control window, the result is invalid. The test should be repeated with a new test card.

Manufacturer:
IMMUNOSPEC CORPORATION
7018 Owensmouth Ave. Suite 103
Canoga Park, CA 91303
Tel : (818)710-1281

PIR13-524
REV.1
Rev date:10/2006

QUALITY CONTROL

1. The control band is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
2. Good Laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control materials which is not provided with this test kit may be commercially available.

LIMITATIONS

1. The test is for qualitative detection of anti-Hantavirus antibody in human serum, plasma or blood sample and dose not indicate the quantity of the antibodies.
2. The test is for *in vitro* diagnostic use only.
3. As in case of all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but should rather be made after all the clinical findings have been evaluated.

PERFORMANCE CHARACTERISTICS

1. Accuracy

Five hundred and sixteen (516) confirmed HFRS samples and four hundred and sixty nine (469) normal samples were tested. The Immunospec Hantavirus Test showed 97.7% accuracy.

<i>Confirmed Clinical Results</i>				
<i>Immunospec Hantavirus Test</i>		Positive	Negative	Total
	Positive	499	6	505
	Negative	17	463	480
Total		516	469	985

Sensitivity = $499/517 = 96.7\%$
 Specificity = $463/480 = 98.7\%$
 Positive predictive value (PPV) = $499/505 = 98.8\%$
 Negative predictive value (NPV) = $463/480 = 96.5\%$
 Accuracy = $(499 + 463)/985 = 97.7\%$

2. Interference

These compounds do not affect the test sensitivity and specificity at the indicated concentrations.

bilirubin 10 mg/dL
 hemoglobin 20mg/dL
 triglycerides 300 mg/dL

REFERENCES

1. Khan A, et al. Viral hemorrhagic fevers, Seminars in Pediatric Infectious Diseases, 8 (suppl.1):64-73, 1997.
2. Linderholm M, Elgh F. Clinical characteristics of hantavirus infections on the Eurasian continent, Current Topics in Microbiology and Immunology, 256: 135-51, 2001
3. Peters CJ, et al. Spectrum of hantavirus infection: hemorrhagic fever with renal syndrome and hantavirus pulmonaary syndrome, Ann. Rev. of Medicine, 50:531-545, 1999
 Schmalijohn C and Hjelle B, Hantaviruses: a global disease problem, Emerging Infectious Diseases, 3:95-104, 1997.
 S:\ImmunoSpec\Immuno doc\Immunospec
 Inserts\Immunospec Rapid test R13-524 Hantavirus



European Authorized Representative:

CEpartner4U , Esdoornlaan 13, 3951DB Maarn
. The Netherlands. Tel.: +31 (0)6.516.536.26

